# Axelgaard

## AG900 SERIES

### AG936 GROUNDING GEL

DESIGNED FOR HIGH FREQUENCY GROUNDING PADS

Property / Test Method	Specification	
Volume Resistivity	7000 ohm-cm max	
Stainless Steel (Skin Side) Adhesion (180° Peel)	110 g/in (43 g/cm) min <sup>1</sup>	
Top Liner (Tie Side) Adhesion to Gel (90° Peel)	200 g/in (79 g/cm) max	
Bottom Liner (Skin Side) Adhesion to Gel (90° Peel)	120 g/in (47 g/cm) max	
Hydrogel Thickness	36 +/- 5 mils (0.9 +/- 0.1 mm)	
Hydrogel pH	4.25 (+/- 0.5)	
Hydrogel Dimensions	Width - 9" (22.9 cm), Length - 300' (91 m) +/-2%	
Top Release Liner	Width - 9.75" (24.8 cm), Thickness - 2 mil (0.05 mm) LDPE	
Bottom Release Liner	Width - 10" (25.4 cm), Thickness - 5 mil (0.13 mm) PET	
Cytotoxicity Test (ISO 10993-5)	Non-Cytotoxic	
Irritation Test (ISO 10993-10)	Non-Irritating	
Sensitization Test (ISO 10993-10)	Non-Sensitizing	
Shelf Life	Slit - 6 months <sup>2</sup> , Non-slit - 12 Months <sup>2</sup>	
Core Size	Slit - 3" (8 cm), Non-slit - 6" (15 cm)	
TYPICAL PROPERTIES - BULK ROLL HYDROGEL		
Property / Test Method	Typical Properties	
Volume Resistivity	4200 ohm-cm	
Stainless Steel (Skin Side) Adhesion (180° Peel)	255 g/in (100 g/cm) <sup>1</sup>	
Top Liner (Tie Side) Adhesion to Gel (90° Peel)	16 g/in (6 g/cm)	
Bottom Liner (Skin Side) Adhesion to Gel (90° Peel)	13 g/in (5 g/cm)	
Hydrogel Thickness	36 mils (0.9 mm)	
Hydrogel pH	4.2	
	0.58	

 $\triangle$  Converter assumes responsibility for material compatibility.

STORAGE CONDITIONS					
Short-Term Temperature (Less than One Month)		32°F - 104°F (0°C - 40°C)			
Long-Term and Optimal Temperature		41°F - 80°F (5°C - 27°C)			
PACKAGING INFORMATION					
Box Size	Box Size		17" x 17" x 12" (43 cm x 43 cm x 30 cm)		
Box Weight		Approx. 68 lbs. (31 kg)			
Pallet Size		Approx. 45" x 45" x 54" (114 cm x 114 cm x 137 cm)			
Pallet Weight (16 Boxes per Pallet)		Approx. 1,118 lbs. (507 kg)			
Headquarters: Fallbrook, CA 92028 USA Europe: DK-8520   Phone: +1 (760) 451-8000 Phone: +45 8674		Lystrup, Denmark 2088	Email: support@axelgaard.com Website: www.axelgaard.com		
Made in the USA	Patented	AG936 12/24		©2024 Axelgaard Manufacturing Co., Ltd.	

# ■Axelgaard

## **AG900 SERIES**

#### **AG936 GROUNDING GEL** DESIGNED FOR HIGH FREQUENCY GROUNDING PADS

### RECOMMENDED STORAGE CONDITIONS FOR BULK ROLL HYDROGELS

Bulk roll hydrogels should be stored at ambient temperature at specific ranges noted below in Number 3a and 3b. Storage conditions outside these ranges can adversely affect product stability and product performance.

#### THE FOLLOWING POINTS SHOULD BE FOLLOWED TO ENSURE ANTICIPATED PRODUCT PERFORMANCE

- 1. Keep hydrogel out of direct sunlight
- 2. Store in the original packaging
- 3. Maintain ambient storage conditions through product expiration
  - a. Short-Term (less than one month) Temperature: 32°F 104°F (0°C 40°C)
    - b. Long-Term and Optimal Temperature: 41°F 80°F (5°C 27°C)
- 4. To prevent dry out of hydrogel, liners should be kept on at all times until the gel is used or converted.

#### SHELF LIFE

The shelf life of bulk roll hydrogels is guaranteed only through the shelf life of the unconverted hydrogel. Axelgaard has stipulated a 12-month shelf life from the date of manufacture for AG936 non-slit bulk roll hydrogel and a 6-month shelf life from the date of manufacture for all slit bulk roll hydrogels.

 $\Delta$  Converting and packaging the AG936 as soon as possible is recommended for optimal processing of the material.

 $\Delta$  Compatibility testing of materials in contact with the hydrogel is also recommended. Certain materials may cause incompatibility when applied to the hydrogel.

#### PACKAGING SPECIFICATION

#### Master Roll Detail

#### Bulk Roll Packaging Slit Gel Roll Packaging 31.5 MSI Units (1 MSI = 1,000 square inches) Plastic Plug **Plastic Plug** ۲ 6" (15.2 cm) 9 3" (7.6 cm) Roll Length: 300' (91 m) +/-2% ò **Diecut Pad Diecut Pad** Roll Width: 10" (25.4 cm) **Release Liner** $^{\circ}$ LDPE Top Liner $\bigcirc$ Slit Gel Roll Width: 9.75" (24.8 cm) 0 Gel Roll **Release Liner** $\circ$ **Double Core Plug** 9 **Diecut Pad** Release Liner $\stackrel{\bullet}{\bigcirc}$ Plastic Plug Ð 0 Slit Gel Roll Max Diameter: 17" (43.2 cm) 6" (15.2 cm) **Release Liner** $^{\circ}$ PET Bottom Liner Tie Layer Shipping Box Width: 10" (25.4 cm) (Wrapped in **Diecut Pad** 0 plastic bag) **Plastic Plug** • 3" (7.6 cm) Skin Layer Shipping Box (Wrapped in Gel Width: 9" (22.9 cm) plastic bag) Headquarters: Fallbrook, CA 92028 USA Europe: DK-8520 Lystrup, Denmark

Phone: +1 (760) 451-8000

Phone: +45 8674 2088

Email: support@axelgaard.com Website: www.axelgaard.com

Patented